


INTRA-AORTIC BALLOON

NEW SIZING GUIDE



Note: This information is to be used as a guidance only. Clinical information and patient factors such as torso length should be considered when selecting the appropriate balloon size. SENSATION and SENSATION PLUS are fiber-optic IAB catheters.

INTRA-AORTIC BALLOON NEW SIZING GUIDE

MAQUET
GETINGE GROUP

Why are the NEW sizes better for my patients?

Larger volume balloons (NEW 30cc, 40cc, 50cc) displace more blood in the aorta during diastole, resulting in improved augmentation and unloading.

With these new sizing guidelines, a larger patient population can benefit from **greater hemodynamic support** than ever before.

Benefits of larger volume IABs

- More blood volume displacement
- More diastolic augmentation
- More systolic unloading

24hr. Clinical Support:

From U.S. and Canada

1 800 777 4222

From outside U.S. and Canada +1 973 244 6100

MAQUET

1300 MacArthur Blvd.

Mahwah, NJ 07430, USA

Phone: +1 201 995 8700

MAQUET Medical Systems USA


45 Barbour Pond Drive

Wayne, NJ 07470, USA

Phone: 1 888 880 2874

<http://ca.maquet.com>

www.maquet.com



Hemodynamic Comparison*

Greater hemodynamic support and higher efficacy for every patient.

* Bench testing completed by MAQUET. Data on file. Bench test results are not necessarily predictive of clinical results.

PN: 0002-08-1092 R3 - (U.S.) MCV00009481 REV.C - (Int'l.) MCV-SB-80000212-EN-02
MAQUET is a trademark of MAQUET GmbH - The following are trademarks ofDatascope
Trademark Corp.: MEGA, SENSATION, LINEAR, SENSATION PLUS. Copyright 2012
Datascope Corp. or its affiliates. All rights reserved. 09/12 - Caution: U.S. Federal Law
restricts this device to sale by or on the order of a physician. ΔRefer to Instructions for
Use for current indications, warnings, contraindications, and precautions.